



SENZA GEN

How are the *in vitro* GARD™ assays used
in real life?

Cases and applications for testing of skin and respiratory sensitizers

Andy Forreryd, PhD

About SenzaGen



- Founded in 2010
- Spin-out from Lund University after over 10 years research
- Highly multidisciplinary team: 20+ employees today



- Business model: global industrial and CRO partnerships
- Own laboratory, continuous development of the technology



Our lead product, GARD™ stands for Genomic Allergen Rapid Detection and is a state-of-the art test platform for assessment of chemical sensitizers



- Launched assays:
- GARD™skin/GARD™potency for skin sensitization
 - GARD™air for respiratory sensitization
 - GARD™skin Medical Device for skin sensitization assessment of medical devices




- Items tested include:
- Active pharmaceutical ingredients
 - Cosmetic ingredients
 - Industrial chemicals
 - Agrochemicals
 - Medical device materials
 - UVCBs, Pre-/pro haptens



- Since Sept. 2017, SenzaGen AB's shares have been traded on Nasdaq First North, Stockholm (SENZA)
- Partners in USA, EU, China, Korea


SENZA
GEN

The **GARD platform** – Currently available assays



GARD™skin (200 genes)

To identify the skin sensitization hazard of chemicals



GARD™skin Medical Device (200 genes)

To identify the skin sensitization hazard of medical devices



GARD™potency (51 genes)


Skin sensitization potency classification according to GHS/CLP



GARD™air (28 genes)

To identify the respiratory sensitization hazard of chemicals

The GARD platform – One testing platform for many answers






The GARD technology platform

Human relevant cells in combination with Genomics and machine learning


INVISIO®

The GARD platform – how it works

SenzaCells: a human dendritic-like cell-line




Cellular responses



The GARD platform – how it works

Cellular response is monitored using biomarker signatures – **Not only a single biomarker**

GARDskin prediction signature 200 genes



Recognition of foreign substances

e.g. **TLRs**, **RXR**, **AHR**



Immunological self-defence mechanisms

e.g. **CD80**, **CD86**



Cellular stress responses


e.g. **NRF2**-pathway



Communication


e.g. chemotaxis receptors

The GARD platform – how it works




How to **GARD** your product - in **6 Steps**


1 GARD Input Finder



Grown SenzaCells




Add different concentrations of the test substance to the cells




Determine the concentration of the test substance where the cells react and 90% survive

2 GARD Main Stimulation



Take test substance at determined concentration and add to fresh batch of cells



Quality control of the cells

3 RNA extraction




Extract RNA from the cells

4 Gene expression profiling



Check the RNA quality




Mix the isolated RNA with reporter probes and load onto a cassette



Analyze the probe cassette to quantify the RNA

5 GARD data analysis application




■ Sensitizer (Train)
■ Non-Sensitizer (Train)
■ Unknown

GARD DV > 0 = Sensitizer
 GARD DV < 0 = Non sensitizer

Upload the results to the GDA web app. One press of the button and the algorithm crunches the data

6 Results



The results are yours!

Dear customer,


We have determined that your test substance is a non-sensitizer.

Kindly,
SENZA GEN


SENZA
GEN

How to **GARD** your product - in **6 Steps**


1 GARD Input Finder



Grown SenzaCells




Add different concentrations of the test substance to the cells




Determine the concentration of the test substance where the cells react and 90% survive

2 GARD Main Stimulation



Take test substance at determined concentration and add to fresh batch of cells




3 RNA extraction



Extract RNA from the cells

4 Gene expression profiling




Check the RNA quality

Mix the isolated RNA with reporter probes and load onto a cassette




Analyze the probe cassette to quantify the RNA

5 GARD data analysis application




Upload the results to the GDAA web app. One press of the button and the algorithm crunches the data

6 Results




How to **GARD** your product - in **6 Steps**


1 GARD Input Finder



Grown SenzaCells




Add different concentrations of the test substance to the cells




Determine the concentration of the test substance where the cells react and 90% survive

2 GARD Main Stimulation



Take test substance at determined concentration and add to fresh batch of cells



Quality control of the cells

3 RNA extraction




Extract RNA from the cells

4 Gene expression profiling



Check the RNA quality




Mix the isolated RNA with reporter probes and load onto a cassette



Analyze the probe cassette to quantify the RNA

5 GARD data analysis application




■ Sensitizer (Train)
■ Non-Sensitizer (Train)
■ Unknown

GARD DV > 0 = Sensitizer
 GARD DV < 0 = Non sensitizer

Upload the results to the GDAA web app. One press of the button and the algorithm crunches the data

6 Results



The results are yours!

Dear customer,


We have determined that your test substance is a non-sensitizer.

Kindly,


SENZA GEN

How to **GARD** your product - in **6 Steps**


1 GARD Input Finder



Grown SenzaCells




Add different concentrations of the test substance to the cells




Determine the concentration of the test substance where the cells react and 90% survive

2 GARD Main Stimulation



Take test substance at determined concentration and add to fresh batch of cells



Quality control of the cells

3 RNA extraction




Extract RNA from the cells

4 Gene expression profiling



Check the RNA quality




Mix the isolated RNA with reporter probes and load onto a cassette



Analyze the probe cassette to quantify the RNA

5 GARD data analysis application




■ Sensitizer (Train)
■ Non-Sensitizer (Train)
■ Unknown

$GARD\ DV > 0 = \text{Sensitizer}$
 $GARD\ DV < 0 = \text{Non sensitizer}$

Upload the results to the GDAA web app. One press of the button and the algorithm crunches the data

6 Results



The results are yours!

Dear customer,

We have determined that your test substance is a non-sensitizer.

Kindly,

SENZA GEN

GARD_{skin} – Performance data

Data set	Sensitivity	Specificity	Accuracy	Reference
Validation of the GARD assay	93% (51/55)	96% (24/25)	94% (75/80)	Johansson, 2019
Accumulated performance	92% (134/145)	81% (50/62)	89% (184/207)	-

TOXICOLOGICAL SCIENCES, 170(2), 2019, 374–381

doi: 10.1093/toxsci/kfz108

Advance Access Publication Date: May 17, 2019

Research Article

Validation of the GARDTMskin Assay for Assessment of Chemical Skin Sensitizers: Ring Trial Results of Predictive Performance and Reproducibility

[Henrik Johansson](#),¹ [Robin Gradin](#),¹ [Angelica Johansson](#),¹ [Els Adriaens](#),² [Amber Edwards](#),³
[Veronika Zuckerstätter](#),⁴ [Anders Jerre](#),¹ [Florence Burleson](#),³ [Helge Gehrke](#),⁴ and [Erwin L Roggen](#)⁵

Within Laboratory Reproducibility:

SenzaGen 82%

Burleson 83%

Eurofins 89%

Between Laboratory Reproducibility:

Concordance: 92%

GARD_{skin} & GARD_{potency} – REACH registration

“The REACH Regulation allows the use of non-adopted in vitro methods in case they meet the EURL ECVAM criteria for entering pre-validation. **For the GARD assay this criteria is met**, as it is currently being validated.”

“The current REACH information requirements require that three KEs are examined and **GARD assay can be used to assess the KE 3.**”


GARD_{skin} – Broad applicability domain

The standard assay protocol is applicable for small molecules in general:

from cosmetics ingredients to various industrial chemicals, including pre- and pro-haptens

SenzaGen and its partners also offer lab services for the explorative assessment of 'difficult-to-test' substances:

- UVCBs
- Natural extracts and mixtures
- Pesticides
- Pharmaceutical ingredients (drug candidates)
- ...and more



Source: SenzaGen in-house studies, excluding data from customer projects

More info at: www.senzagen.com

[Home](#)[Services](#)[Science](#)[About](#)[Partners](#)[Investors](#)[News and Events](#)[Contact](#)[GARD® Technology](#)[Science Behind](#)[Scientific Publications](#)[Posters](#)[Application Notes](#)[Webinars](#)

NEXT GENERATION SAFETY TESTING.

SenzaGen's tests identify potentially harmful chemicals. It is more relevant and cost-efficient, making animal testing unnecessary and daily products safer.



NEW PRODUCT LAUNCH

GARD™skin Medical Device

At Eurotox Sep 8-11, Helsinki



Case study I – GARD_{skin} for "difficult to test" substances

Why substances can be difficult to test:

Pre/Pro haptens: Activation may be needed to create the allergen.

SenzaCells:

Aldehyde dehydrogenases (ALDH)

Cytochrome p-450 (CYP)

N-acetyltransferase 1 (NAT-1)

Solubility: High log $P_{o/w}$ value or other properties leading to low solubility in aqueous media.

Case study I – GARD_{skin} for "difficult to test" substances

Compound	Pre/pro hapten	High logP	DPRA (TG442C)	ARE-NRF2 (TG 442D)	h-CLAT (TG442E)	GARD (TGP 4.106)
2-Aminophenol	YES	-	S	S	S	S ¹
2-nitro-1,4-Phenylendiamine	YES	-	S	S	S	S ¹
Aniline	YES	-	NS	NS	S	S ¹
Cinnamic alcohol	YES	-	S	S	S	S ¹
Ethylene diamine	YES	-	NS	S	S	NS ¹
Eugenol	YES	-	S	NS	S	S ¹
Geraniol	YES	-	NS	S	S	S ¹
Hydroquinone	YES	-	S	S	S	S ²
Isoeugenol	YES	-	S	S	NS	S ¹
Linalool	YES	-	NS	NS	S	S ¹
p-Phenylenediamine	YES	-	S	S	S	S ¹
Propyl gallate	YES	-	S	S	S	S ¹
Resorcinol	YES	-	NS	NS	S	S ¹
Farnesol	YES	-	NS	S	S	S ¹
Abietic acid	YES	YES (6.5)	S	S	NS	S ¹
Chlorpromazine	YES	YES (5.4)	NA	NS	S	S ¹
Lauryl gallate	YES	YES (6.2)	S	S	S	S ³
Amylcinnamyl alcohol	YES	YES (4.4)	S	NS	NS	S ¹
Limonene	YES	YES (4.6)	NS	NS	S	S ¹
Benzoyl peroxide	-	YES (3.5)	NS	NS	S	S ³
Hexylcinnamic aldehyde	-	YES (4.8)	S	NS	NS	S ¹
Isopropyl myristate	-	YES (8.5)	NS	NS	S	NS ⁴
propyl paraben ^{*1}	-	YES (3.4)	NS	S	S	S ¹
Tocopherol	-	YES (6.9)	NS	S	NS	S ¹
Accuracy			61%	58%	71%	92%

References:

¹ Johansson et al. 2017

² Forreryd et al. 2016

³ Zeller et al. 2017

⁴ Johansson et al. 2019

*Basketter

Human potency Class 5


SENZA
GEN

Case study II – Testing of UVCBs


UVCBs: Unknown or Variable composition, Complex reaction products and Biological materials

The challenges: Complexity and uncertainty of compositions; Hydrophobic, often has solubility issues for water-based test systems





Vehicle: DMSO



Vehicle: H₂O



Vehicle: DMF/Glycerol



Available solvent: Acetone, Chlorobenzene, DMF, DMSO, Isopropanol, Ethanol, Glycerol, Hexane, Olive oil and sesame oil.

Case study II – Testing of UVCBs

Materials & Methods

Test item	GARDskin
7 UVCBs with known average MW	Additional solubility tests performed to select suitable solvents

Results in comparison with existing *in vivo* data

Test items	Existing data	Existing classification	GARDskin
1	LLNA, Buehler	S, 1B	S
2	LLNA	NS	S
3	LLNA, Buehler	S, 1B	S
4	LLNA	S, 1B	S
5	Buehler, HRIPT	NS	NS
6	LLNA, Buehler, HRIPT	S, 1B	S
7	Buehler	S, 1B	S

Conclusion

GARDskin:

good applicability potential for UVCBs, consistent with *in vivo* data.

The assessment of UVCBs need to be handled **case by case**.


Case study III– Testing of natural mixtures

Henna: natural dye from *Lawsonia inermis*.

Henna-based hair colouring products: often considered safer alternatives to synthetic hair dyes containing known skin sensitizers such as p-phenylenediamine (PPD). However, skin allergic reactions to henna products have been reported.



Materials & Methods



*The images are from the internet and only serve as examples of commercially available henna products

This study is a collaboration between Federal University of Goiás, Brazil and Lund University, Sweden. For more details, please see the original publication (de Ávila et al. Contact Dermatitis, 2019)

Case study III– Testing of natural mixtures

Table 1. GARDskin prediction results for commonly used hair dye ingredients

Test materials	Human classification	Animal classification	mDPRA	IL-18	USENSTM	GARDTM skin
Reference controls						
Dimethyl sulfoxide	NS	NS	NS	NS	NS	NS
Glycerol	NS	NS	NS	NS	NS	NS
Sodium dodecyl sulfate	NS	S	S	NS	S	NS
DNCB	S	S	S	S	S	S
Eugenol	S	S	S	S	S	S
2-Hydroxyethyl acrylate	S	S	-	-	-	S
Hair dye ingredients						
1,4-Diaminoanthraquinone	S	S	S	S	S	S
2-Amino-3-hydroxypyridine	S	NS	S	S	S	S
Lawson	NA	Equivocal	S	NS	S	S
5-Amino- <i>o</i> -cresol	S	S	S	NS	S	S
Hydroquinone	S	S	S	S	S	S
<i>p</i> -Phenylenediamine	S	S	S	S	S	S
Resorcinol	S	S	S	S	S	S
Disperse orange 3	S	S	S	S	S	S
Basic red 51	NA	NS	S	S	S	S
Pyrogallol	S	NS	S	S	S	S
Concordance vs. Human data			91.7%	92.3%	92.3%	100%
Concordance vs. Animal data		78.5%	76.9%	64.3%	78.6%	73.3%

Case study III– Testing of natural mixtures

Table 2. GARDskin prediction results for 10 commercial henna-based products

Henna Product	PPD (label)	PPD	mDPRA	IL-18	USENS™	GARDskin™
1		1.091 ± 0.028	S	S	S	S
2	+	2.970 ± 0.046	S	S	S	S
3		0.030 ± 0.001	S	S	NS	S
4		0.032 ± 0.006	S	S	NS	S
5		4.321 ± 0.028	S	S	S	S
6		1.020 ± 0.100	S	S	S	S
7		0.577 ± 0.015	S	NS	S	S
8	+	2.541 ± 0.057	S	NS	S	S
9		0.760 ± 0.017	S	S	S	S
10		3.354 ± 0.163	S	S	S	S

Conclusion

GARDskin show high concordance to human data for testing of **hair dye ingredients**

GARDskin is applicable for **henna-based hair colouring products (powder mixtures)**.

Case study IV – Testing of oil-based mixtures

The Customer: a leading supplier of natural based functional ingredients/ raw materials for cosmetic and personal care applications


The problem: An oil-based mixture with **unpleasant and characteristic smell**

- Used as raw materials sold to cosmetic and personal care companies
- The oil-based mixture contains **chemical x** as contaminant, which is the suspect for the unpleasant smell.
- Distillation to get rid of contaminants of **chemical x**

Case study IV – Testing of oil-based mixtures

Methods & Results (GARDskin)

Test item 1	Test item 2	Test item 3
Pure chemical x	Oil mixture containing 1200 ppm chemical x	Distilled oil mixture containing 10 ppm chemical x
Non sensitizer	Non sensitizer	Sensitizer



Conclusion

This case is a good example of how GARDskin can be used for assessment of oil-based mixtures and essential oils.



Case study V– Testing Medical Devices for skin sensitization hazard

***In vitro* skin sensitization testing in the Biological Evaluation of Medical Devices according to ISO 10993-10:2010**

- Use polar and non-polar extraction vehicles according to ISO 10993-12:2012 (saline, sesame oil, super refined olive oil and cell culture media)
- Classify leachables in extracts of medical devices as either skin sensitizers or non-sensitizers
- 200 genomic biomarkers in GARDskin Predictive Signature (GPS)

Case study V– Testing Medical Devices for skin sensitization hazard

Extraction steps of solid material follow the international standards for safety assessment of **Medical Device (ISO 10993-12)**.



Case study V– Testing Medical Devices for skin sensitization hazard

Table 1. Summary of the results from the in-house validation of GARD®skin Medical Device compared with LLNA (as listed in the CE STTF database) and Human potency classification (HP) for the chemicals (Basketter et al. 2014) .

Test material	Chemical	Sensitizing potential		GARD®skin Medical Device Prediction		
		LLNA	HP	Saline	Olive oil	Sesame oil
Silicone	None	N/A	N/A	Non-sensitizer	Non-sensitizer	Non-sensitizer
	2-aminophenol	Strong	Cat 2	Sensitizer	Sensitizer	Sensitizer
	Cinnamic aldehyde	Moderate	Cat 2	Sensitizer	Sensitizer	Sensitizer
	Propyl gallate	Strong	Cat 2	Sensitizer	Sensitizer	Sensitizer
	Phenyl benzoate	Weak	Cat 3	Sensitizer	Sensitizer	Sensitizer
TPU	None	N/A	N/A	Non-sensitizer	Non-sensitizer	Not tested
	Propyl gallate	Strong	Cat 2	Sensitizer	Sensitizer	Not tested
	Phenyl benzoate	Weak	Cat 3	Sensitizer	Sensitizer	Not tested
Silicone tube	-	N/A	N/A	Non-sensitizer	Non-sensitizer	Non-sensitizer
TPU tube	-	N/A	N/A	Non-sensitizer	Non-sensitizer	Non-sensitizer
PVC tube	-	N/A	N/A	Non-sensitizer	Non-sensitizer	Non-sensitizer
Vehicle control	-	Neg		Non-sensitizer	Non-sensitizer	Non-sensitizer
Positive control	p-Phenylenediamine	Pos		Sensitizer	Sensitizer	Sensitizer

Case study VI – Safety assessment of e-liquids

Materials & Methods

In this customer study, **GARDskin**, **GARDpotency** and **GARDair** were used to predict and compare the skin and respiratory sensitisation potential of **three experimental** and **two commercial e-liquids**. To our knowledge this is the first published study assessing e-liquids using in vitro assays.

Regulatory Toxicology and Pharmacology 103 (2019) 158–165



Contents lists available at ScienceDirect

Regulatory Toxicology and Pharmacology

journal homepage: www.elsevier.com/locate/yrtph



The use of Genomic Allergen Rapid Detection (GARD) assays to predict the respiratory and skin sensitising potential of e-liquids




Matthew Stevenson^{a,*}, Lukasz Czekala^a, Liam Simms^a, Nicole Tschierske^a, Olivia Larne^b, Tanvir Walele^a

^a Imperial Brands PLC, 121 Winterstoke Road, Bristol, BS3 2LL, UK

^b SenzaGen, Scheelevägen 2, 22381, Lund, Sweden


Case study VI – Safety assessment of e-liquids

The typical composition of e-liquids



E-liquids are usually **complex mixtures**, especially when they include natural flavourings

Case study VI – Safety assessment of e-liquids





GARDair prediction

A controls: 5 out of 7 accurately classified as respiratory sensitizers, no false positives. Sensitivity 71%, Specificity 100%, Overall accuracy 89%

B test substances: None of the base liquids nor commercial e-liquids were classified as respiratory sensitizers.

Case study VI – Safety assessment of e-liquids



GARDskin prediction :

2 commercially available flavoured e-liquids were predicted as skin sensitizers

GARDpotency prediction :

Test substances: 2 commercially available flavoured e-liquids were further classified as weak skin sensitizers (1B)

Case study VI – Safety assessment of e-liquids

Composition of test materials. PG: Propylene glycol; VG: Vegetable glycerol; BL: base liquid; CF: commercially available flavoured e-liquid.

Test material	Content [w/w %]			
	PG	VG	Nicotine	Other substances (e.g. water, flavourings)
BL 0% ^a	50	50	–	–
BL 1.6% ^a	49.2	49.2	1.6	–
BL 4.5% ^a	47.75	47.75	4.5	–
CF Blu Cherry 1.6% ^a	40.4	53.4	1.6	4.6
CF 1.2% ^a	66.5	18.8	1.1	13.6

^a % refer to % nicotine content.

Test substance classification with the GARD assay for respiratory sensitisers, GARDskin and GARDpotency assays.

	Respiratory assay	GARDskin	GARDpotency
BL 0%	Non-Sensitiser	Non-Sensitiser	No Cat
BL 1.6%	Non-Sensitiser	Non-Sensitiser	No Cat
BL 4.5%	Non-Sensitiser	Non-Sensitiser	No Cat
CF Blu Cherry 1.6%	Non-Sensitiser	Sensitiser	1 B
CF 1.2%	Non-Sensitiser	Sensitiser	1 B

Conclusion

The **GARD platform** offers an integrated test strategy for assessment of skin and respiratory sensitization potential of **complex mixtures** such as e-liquids.

The assays are potentially useful for product development and hazard identification of other types of complex formulations.

The **GARD platform** – How can I get my **substances tested?**

Contact us:

We help to design the testing strategy for your specific substances.

Select assay(s):

GARDskin, GARDpotency, GARDair, GARDskin Medical Device

Test Substance Questionnaire:

Do you have information on vehicles?
If not, we evaluate it for you.

Turnaround time:


4-6 weeks.

Sample requirements:


0.5 g (solids) or 1 ml (liquids). Can be adapted to lower amounts.

Your preferred testing site:

Select where you like to have your testing done.



The **GARD platform** – Where can I get my **testing done?**



• **Licence Labs*:**

Burleson Research Technologies
Eurofins BPT
MB Research Laboratories

• **Distributors*:**

Charles River Laboratories
Eurosafte
Guangzhou CHN-ALT Biotech Co., Ltd
PKDerm
Woo Jung BSC
XCellR8

* In alphabetic order

Information and contacts

www.senzagen.com

SENZA GEN

[Home](#) [Services](#) [Science](#) [About](#) [Partners](#) [Investors](#) [News & Events](#) [Contact](#)

ANIMAL FREE TESTING.

Our method mimics the human immune system.

[OUR SERVICES](#) [LEARN MORE](#)

SenzaGen rings the bell
to open the stock market and the start for trading in #SENZA at Nasdaq First North.

"We are excited to join Nasdaq First North and are proud of the support from investors whom have trusted us and joined our journey as a listed company," said Anki Malmberg Hager, CEO of SenzaGen.

"We welcome SenzaGen to our European growth market, Nasdaq First North," said Adam Kostyál, SVP and Head of European listings at Nasdaq. "SenzaGen operates in an exciting field and we look forward to supporting them for many years to come."

Andy Forreryd, PhD

Key Account Manager & Scientific Liaison

andy.forreryd@senzagen.com

+46 734331777

Joshua Schmidt, PhD

Business Development Director Americas

Joshua.schmidt@senzagen.com

+1 651.440.5691

SENZA GEN