

## GARD<sup>®</sup>skin Dose-Response

*In vitro* quantitative potency assessment



**GARDskin Dose-Response is an *in vitro* test for quantitative skin sensitizing potency assessment of chemicals, adapted from GARDskin (OECD TG 442E).**

GARDskin Dose-Response provides an estimated threshold concentration for a test substance to induce skin sensitizing effects; a lower concentration equals a higher potency and vice versa.

This threshold concentration can be used for skin sensitizing potency ranking of candidate ingredients and early decision-making during product development.

Additionally, readouts provide a prediction of correlating LLNA EC3 values, which is traditionally used to measure skin sensitizing potency of chemicals. Furthermore, the readout can also predict human potency and/or GHS/CLP classification, all with high statistical significance.

**Non-animal readout to predict:**

- LLNA EC3 value
- Human potency
- GHS/CLP classification

### Features and Benefits

#### Test system

- Human dendritic-like cell line: SenzaCell™.

#### Solvent

- Standard: DMSO and H<sub>2</sub>O.
- Other available solvents: Acetone, DMF, Isopropanol, Ethanol, Glycerol, Olive oil, Sesame oil.

#### What it measures

- GARDskin response values in a titrated range of multiple concentrations.

#### Readout

- Hazard: skin sensitizer or non-sensitizer.
- Quantitative potency: cDV<sub>0</sub>.\*
- Further prediction:
  - LLNA EC3 value.
  - Human potency.
  - GHS/CLP classification 1A or 1B.

#### High performance

- The performance of GARDskin Dose-Response is associated with GARDskin, which demonstrated a 94% accuracy for skin sensitizing hazard prediction.\*\*
- GARDskin Dose-Response readouts have shown significant correlation with the LLNA EC3 values and human potency.

#### Short turnaround time

- 6-8 weeks for standard studies.

#### Low sample requirement

- 0.5g (solid) or 1ml (liquid).

#### Compliance

- GARDskin Dose-Response is an adaptation of GARDskin, which is approved by OECD as part of Test Guideline 442E for *in vitro* skin sensitization.
- GLP or non-GLP.

\*cDV<sub>0</sub> is the threshold concentration for the test substance to induce a positive response in GARDskin.

\*\*Johansson et al., Validation of the GARDskin assay for assessment of chemical skin sensitizers – ring trial results of predictive performance and reproducibility. Toxicological Sciences. May 17, 2019.

## How GARD® works

The GARDskin Dose-Response investigates the GARDskin response values in a titrated range of multiple concentrations in a dose-response manner to find the threshold concentration required to induce a positive decision value.

This threshold concentration,  $cDV_0$ , is used to estimate the inherent sensitizing potency of a chemical. A low value indicates a high inherent skin sensitizing potency and a high value indicates a low inherent skin sensitizing potency.

$cDV_0$  can also be used to predict the traditional EC3 value from Local Lymph Node Assay, the human skin sensitizing potency, and/or GHS/CLP classification 1A or 1B depending on the need.

All the GARD assays are based on the same technology platform. Read more about the GARD technology platform and assay development principles on [www.senzagen.com/science](http://www.senzagen.com/science).


## In-house validation study

### Method

An in-house validation study of GARDskin Dose-Response was conducted to investigate the dose-response relationship between GARDskin classifications and test chemical concentration. Data was generated on 29 reference chemicals by performing the GARDskin assay in a titrated range of concentrations.

### Results

The results confirmed that  $cDV_0$  informs on the skin sensitizing potency and can be used to directly monitor sensitizing potency of chemicals.  $cDV_0$  values from skin sensitizers were associated with GHS classification labels, as well as significantly correlated to both human and LLNA potency reference data (Figure 2).



	GARD	LLNA
Response value	DV	SI
Binary Threshold	DV=D	SI=3
Readout	$cDV_0$ (µM)	EC3 (%)

Figure 1. The experimental setup of GARDskin Dose-Response allowing for establishment of the  $cDV_0$  value, which is derived analogously to the LLNA EC3 value.




Figure 2. Scatter plots displaying the relationship between estimated  $cDV_0$  values and A) LLNA EC3 values and B) human NOEL values. The fitted lines represent linear regression models fitted to the data, and the shaded areas describe the 95% confidence intervals of the fits. Encircled datapoints indicate pre- and pro-haptens.

